



Muskegon Laboratory: 231-672-6259

WHAT COLOR TUBE TO DRAW

ALL TEST LISTED CAN BE TRANSPORTED TO LAB VIA ROUTINE COURIER	
Tubes with gel in them should be spun within 2 hours and prior to transport	
ALBUMIN LAB45	INTERFERON GAMMA FOR TB LAB4702 -NEED KIT
ALKALINE PHOSPHATASE LAB112	IONIZED CALCIUM -SHORT LITH GREEN LAB54
ALT LAB132	IRON LAB94
AST LAB131	IRON & IRON BINDING STUDIES, TIBC, %SAT LAB829
AMYLASE LAB48	LEAD LAB98 **FORM NEEDED**
ANTINUCLEAR ANTIBODY (ANA) LAB4147	LUTENIZING HORMONE LAB87
ANTINUCLEAR ANTIBODY (ANA) REFLEX LAB4146	LIPASE LAB99
BASIC METABOLIC PANEL (BMP) LAB15	LIPID PANEL LAB18
BETA (QUANT) HCG LAB143	LITHIUM LAB29
BILIRUBIN DIRECT LAB52	MAGNESIUM LAB103
BILIRUBIN TOTAL LAB50	METHYLMALONIC ACID LAB835
BLOOD TYPE & ANTIBODY SCREEN -TALL LAV LAB276	MONO LAB482
BLOOD UREA NITROGEN (BUN) LAB140	MUMPS IgG LAB160
B-TYPE NATRIURETIC PEPTIDE (BNP)- TALL LAV LAB106	NOCOTINE AND METABOLITES- TALL RED LAB4626
CALCIUM LAB53	PARATHYROID, INTACT- TALL LAV LAB108
CARBAMAZEPINE (TEGRETOL) LAB21	PHENYTOIN (DILANTIN) LAB31
CARBON DIOXIDE LAB55	PHENYTOIN (DILANTIN) FREE LAB32
CARCINOEMBRYONIC AG LAB57	PHOSPHURUS LAB113
CBC LAB294	POTASSIUM LAB114
CBC WITH DIFF LAB293	PROLACTIN LAB531
CELIAC DISEASE PANEL LAB4254	PROTEIN LAB118
CHOLESTEROL LAB860	PSA FREE & TOTAL LAB171
CHLORIDE LAB59	PSA SCREENING LAB116
CLOZAPINE LAB4308	PT-INR LAB320
COMPREHENSIVE METABOLIC PANEL (CMP) LAB17	PTT LAB325 Deliver to lab within 4 hours of draw
COPPER LAB817 **NOT PIPETTE** DK BLUE PLAIN	RENAL FUNCTION PANEL LAB19
CORTISOL LAB61	RETICULOCYTE COUNT LAB296
C-REACTIVE PROTEIN (CRP) LAB149	RHEUMATOID FACTOR LAB206
C-REACTIVE PROTEIN (CRP) HIGH SENSITIVITY LAB150	RUBELLA IgG LAB496
CREATININE LAB383	RUBEOLLA IgG LAB657
CREATINE KINASE LAB62	SED RATE LAB322 SHORT LAV
CYCLOSPORINE LAB874	PROTEIN ELECTROPHORESIS SERUM LAB119
CYSTATIN C LAB4343	SIROLIMUS (RAPAMUNE) LAB875
D-DIMER LAB313 Deliver to lab within 4 hours of draw	SODIUM LAB122
DHEA LAB522	SYPHILIS ANTIBODY LAB8418
DHEA SULFATE LAB524	TACROLIMUS (PROGRAF) LAB876
DIGOXIN LAB23	TESTOSTERONE LAB124
ELECTROLYTES LAB16	TESTOSTERONE FREE LAB5888
FERRITIN LAB68	TESTOSTERONE FREE, BIOAVAILABLE AND TOTAL LAB4771-TALL RED
FOLATE LAB69	THYROGLOBULIN ANTIBODY LAB515
FREE T3 LAB137	THYROGLOBULIN LEVEL LAB533
FREE T4 LAB127	THYROID PEROXIDASE ANTIBODY LAB858
FOLLICLE STIMULATING HORMONE LAB86	T3 TOTAL LAB136
GAMMA GLUTAMYLTRANSFERASE (GGT) LAB85	T4 TOTAL LAB126
GENTAMYCIN LAB27	TROPONIN I LAB7518
GLUCOSE LAB82	TSH LAB129
HCG (QUAL) LAB144	TSH REFLEX TO FT4 & FT3 LAB6894
HEMOCHROMATOSIS LAB833	URIC ACID LAB141
HEMOGLOBIN A1C LAB90	VALPROIC ACID (DEPAKOTE) LAB24
HEPATIC FUNCTION PANEL LAB20	VALPROIC ACID (DEPAKOTE) FREE- TALL RED LAB25
HEPATITIS B ANTIBODY LAB472	VANCOMYCIN LAB40
HEPATITIS B ANTIGEN LAB471	VARICELLA IgG LAB162
HEPATITIS PANEL, ACUTE LAB551 SST&LAV	VITAMIN B12 LAB67
HEPATITIS C AB SCREENING LAB4450 SST&LAV	VITAMIN D 1, 25-DIHYDROXY LAB536
HIV SCREEN LAB4477 SST/LAV	VITAMIN D 25-HYDROXY LAB535
IMMUNOFIXATION LAB174	ZINC LAB581 DK BLUE EDTA

DRAW TUBES IN THE FOLLOWING ORDER	
1. BLOOD CULTURE	7. LITH HEP (TALL LIGHT GRN)
2. BLUE TOP	8. SHORT AND TALL LAV
3. RED TOP	9. PINK TUBE
4. PLAIN DK BLUE	10. EDTA DK BLUE
5. SST/GOLD TOP	11. ACD-A OR ACD-B (FAT YELLOW)
6. SODIUM HEPARIN (FAT TALL GRN)	12. QUANTIFERON TB KIT